Acetabulumfräser
Acetabular Reamers

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Directions for use and safety instructions for Acetabular Reamers

Acetabular Reamers may only be implemented by personnel who by virtue of their training and experience are familiar with the use of the said instruments.

Scope of application

The purpose of Acetabular Reamers is to shape the bone to insertion of hip cups. The operator is responsible for selection of the appropriate instrument on each occasion.

Appropriate usage

- Use instruments only in accordance with this information sheet and instructions of the motor manufacturer
- Use only drives which are technically and hygienically perfect, maintained and cleaned
- Use only instruments which are in perfect technical and sterile condition
- Before commissioning check safe fit position ► danger of injury
- When changing tools, secure the drive system in order to avoid accidental activation by foot pedal etc. ► danger of injury
- Bring revolutions up to speed before applying to object – tilting and levering is to be avoided ► danger of breaking instrument
- Damaged instruments (e.g. broken or blunt cutter edges; cracks; corrosion; buckling) may not be used ► overheating and damage to tissue and bone (necrosis)
- Wear appropriate protective clothing (e.g. goggles) during use
- Contact with the tool's cutting edge is to be avoided ► Danger of contamination and infection where operation gloves are damaged
- Contact with the cover is to be avoided ► danger of injury; damage to instrument
- Instruments not to be used for working on metals (e.g. implants) ► danger of injury and infection through abrasion

Gripping pressure

- Excessive grip pressure is to be avoided ► damage to the part (in extreme cases: breaking of tool), excessive heat generation

Cooling

- Take care to cool sufficiently (e.g. by using physiological saline solution) ► insufficient cooling results in a high heat generation, which can in some cases lead to irreversible bone and/or tissue damage (necrosis)

Cleaning; disinfection; sterilisation; storage

- **The instruments are supplied in a non-sterile state**
- Sterilise with damp heat (steam) in an autoclave at 134° C, leaving for at least 18 minutes (please refer to the manufacturer's instructions)
- Temperatures of over 180° C are to be avoided ► loss of the part's hardness = reduction of durability
- All packaging containing sterilised goods should have a sticker attached showing date of sterilisation and expiry date
- Prior to use for the first time, storage should be in the original packaging at room temperature and protected from dust and damp
- Prior to use for the first time and immediately after each use, the instruments must be disinfected, cleaned and sterilised ► danger of infection
- Wear gloves when working with unclean instruments ► danger of infection
- When the instruments are very dirty, cleaning should be carried out in an ultrasound bath with a suitable cleaning substance (see the manufacturer's instructions)
- The instruments should not come into contact with one another during cleaning/disinfection ► risk of damage
- Immediately remove traces of blood, secretion, tissue and bone – dried traces cause corrosion
- Use disinfectant with corrosion protection
- Immediately remove traces of disinfectant and cleaning agents with running water and dry carefully ► danger of corrosion
- Thermal disinfection is not suitable for rotating instruments
- Contact with hydrogen peroxide (H₂O₂) is to be avoided
- Immediately prior to use, instruments should be kept in appropriate, hygienically maintained containers/trays
- During storage, instruments should be protected against damp, dust and contamination

Further information is available from the "Arbeitskreis Instrumentenaufbereitung (AKI)" (Instrument Preparation Working Group) in its brochure "Instrumentenaufbereitung richtig gemacht" (Instrument Preparation done properly) (www.a-k-i.org)

Safety and Liability

- Before use, the user is obliged to check as he sees fit the suitability and application possibilities for the intended purpose – the use of the instrument is the responsibility of the user

RISA GmbH accepts no liability for damage caused by use of the product