



Sinterdiamanten
Sinterdiamonds

Anwendungsbereich

Zum Bearbeiten von Keramik, Chrom-Kobalt, Chrom-Nickel, Aufbrennlegierungen und Gold.

- Vorteile:**
- lange Standzeit und hohe Kantenstabilität.
 - gleichmäßige Schleifergebnisse, bei geringer Wärmeentwicklung.
 - vibrationsfreies Arbeiten.

Wichtiger Arbeitshinweis:

Neue Instrumente, und Instrumente deren Schleifleistung nachlässt sowie Sinterdiamanten die auf der Keramik schwärzen, müssen unbedingt für circa 15 Sekunden bei 1.000 bis 5.000 U / min. an dem mit Wasser angefeuchteten Reinigungsstein gesäubert werden.

Alle Diamantinstrumente sind als HP oder RA (Winkelstück) erhältlich.

Field of application:

For machining of ceramic, chrome-cobalt, chrome-nickel, burnt-on alloys and gold

- Advantages:**
- longevity - ultimate stability of the cutting edges.
 - excellent grinding results - much lower heat generation.
 - vibration-free working.





Important Note:





New instruments, and instruments with reduced grinding power, or Sinterdiamonds blackening the ceramic, have to turn for 15 sec. at 1.000 to 5.000 Rpm at the damped dressing stone.

All instruments are available in HP or RA (contra angle) .

Körnungen

Grit

-  roter Ring = fein
-  ohne Ring = mittel
-  grüner Ring = grob
-  schwarzer Ring = extra grob

-  red ring = fine
-  without ring = medium
-  green ring = coarse
-  black ring = extra coarse

Bestellbeispiel

Ordering example:

Fein / fine: M 801/029F HP
Mittel / medium: M 801/029 HP

Grob / coarse: M 801/029G HP
Extra grob / extra coarse: M 801/029XG HP

Verpackungseinheit / Unit: 1 Stück / pieces / Packung / packet.



Sinterdiamanten
Sinterdiamonds

M 801 021	M 801 029	M 801 035	M 804 016	M 804 018	M 804 027	M 804 037	M 804 047	M 804 080	M 804S 018	M 804S 027
M 804S 037	M 804S 047	M 804S 080	M 805 016	M 805 018	M 805 027	M 805 029	M 805 037	M 805 047	M 805 060	M 805 080
M 805S 016	M 805S 018	M 805S 027	M 805S 037	M 805S 047	M 805S 080	M 807 029	M 807 050	M 809 047	M 809 080	M 809 160
M 809 220	M 818 080	M 825 047	M 825 080	M 830 016	M 835 010	M 837 023	M 837 027	M 837 029	M 837 031	M 837 050



Sinterdiamanten
Sinterdiamonds



M 837	M 837L	M 837L	M 838	M 838	M 838	M 838	M 838	M 838	M 838	M 838	M 838
060	023	031	016	021	023	029	031	050	060	080	080



M 838L	M 838L	M 845	M 845	M 845	M 845L	M 845L	M 846	M 846	M 846L	M 856
023	031	016	031	050	031	050	027	031	031	021



M 856	M 856	M 856L	M 856L	M 860	M 860	M 862	M 862	M 862	M 862	M 863
031	040	023	031	016	031	029	031	050	080	050



M 871	M 879	M 879
050	050	060



Reinigungsstein / Dressing Stone
Art.-Nr.: M 750