



Diamant-Scheiben galvanisch gebunden
Diamond Discs galvanic bound



914 einseitig
914 single-sided



915 doppelseitig
915 double-sided



916 einseitig
916 single-sided



916 einseitig
916 single-sided



917 doppelseitig
917 double-sided





918 doppelseitig
918 double-sided

Beschreibung: Diamant-Scheiben galvanisch gebunden.
Anwendungsbereich: Zum Trennen, Separieren und seitlichen Beschleifen von Keramik, Kunststoffen, Composite, Gips, und Metallen.

Description: Diamond Discs galvanic bound.
Field of application: For cutting, separating, and lateral grinding of ceramic, plastic, composite, plaster and metals.

Bestellnummer / Order Number			Ø 19 mm	Ø 22 mm
Montiert mounted	unmontiert unmounted			
HP 914/190 F	914/190 F	(fein, einseitiger Belag; fine, single-sided coating)	<input type="checkbox"/>	
HP 914/220 F	914/220 F	(fein, einseitiger Belag; fine, single-sided coating)		<input type="checkbox"/>
HP 914/190 XF	914/190 XF	(x-fein, einseitiger Belag; x-fine, single-sided coating)	<input type="checkbox"/>	
HP 914/220 XF	914/220 XF	(x-fein, einseitiger Belag; x-fine, single-sided coating)		<input type="checkbox"/>
HP 915/190 F	915/190 F	(fein, doppelseitiger Belag; fine, double-sided coating)	<input type="checkbox"/>	
H 915/220 F	915/220 F	(fein, doppelseitiger Belag; fine, double-sided coating)		<input type="checkbox"/>
HP 915/190 XF	915/190 XF	(x-fein, doppelseitiger Belag; x-fine, double-sided coating)	<input type="checkbox"/>	
HP 915/220 XF	915/220 XF	(x-fein, doppelseitiger Belag; x-fine, double-sided coating)		<input type="checkbox"/>
HP 916/190 F	916/190 F	(fein, einseitiger Belag; fine, single-sided coating)	<input type="checkbox"/>	
HP 916/220 F	916/220 F	(fein, einseitiger Belag; fine, single-sided coating)		<input type="checkbox"/>
HP 916/190 M	916/190 M	(mittel, einseitiger Belag; medium, single-sided coating)	<input type="checkbox"/>	
HP 916/220 M	916/220 M	(mittel, einseitiger Belag; medium, single-sided coating)		<input type="checkbox"/>
HP 917/190 F	917/190 F	(fein, doppelseitiger Belag; fine, double-sided coating)	<input type="checkbox"/>	
HP 917/220 F	917/220 F	(fein, doppelseitiger Belag; fine, double-sided coating)		<input type="checkbox"/>
HP 917/190 M	917/190 M	(mittel, doppelseitiger Belag; medium, double-sided coating)	<input type="checkbox"/>	
HP 917/220 M	917/220 M	(mittel, doppelseitiger Belag; medium, double-sided coating)		<input type="checkbox"/>
HP 918/220 F	918/220 F	(fein, doppelseitiger Belag; fine, double-sided coating)		<input type="checkbox"/>





Sinterdiamant-Scheiben / Nonflex
Sinterdiamond Discs / Nonflex

Beschreibung: Sinterdiamant-Scheiben / Nonflex

Description: Sinterdiamond discs / Nonflex



Bestellnummer Order Number	Durchmesser Ø Diameter Ø	Stärke Strength	Anwendungsbereich Field of application
M 941/200 HP (mittel, medium)	Ø 20 mm	0,27 mm	Keramik; ceramic
M 941G /220 HP (grob, coarse)	Ø 22 mm	0,30 mm	Metall, Keramik; metal, ceramic
M 941XG /220 HP (x-grob, x-coarse)	Ø 22 mm	0,30 mm	Metall, Keramik; metal, ceramic

Sinterdiamant-Scheiben / Mediumflex
Sinterdiamond Discs / Mediumflex

Beschreibung: Sinterdiamant-Scheiben / Mediumflex

Description: Sinterdiamond discs / Mediumflex




Bestellnummer Order Number	Durchmesser Ø Diameter Ø	Stärke Strength	Anwendungsbereich Field of application
M 942a /210 HP (mittel, medium)	Ø 21 mm	0,20 mm	Kunststoff, Gold; plastic, gold
M 942b /210 HP (mittel, medium)	Ø 21 mm	0,23 mm	Kunststoff, Gold; plastic, gold
M 942 /300 HP (mittel, medium)	Ø 30 mm	0,30 mm	Gips; plaster
M 942 /400 HP (mittel, medium)	Ø 40 mm	0,30 mm	Gips; plaster
M 942 /450 HP (mittel, medium)	Ø 45 mm	0,30 mm	Gips; plaster

Sinterdiamant-Scheiben / Flex
Sinterdiamond-Discs / Flex

Beschreibung: Sinterdiamant-Scheiben / Flex

Description: Sinterdiamond discs / Flex




Bestellnummer Order Number	Durchmesser Ø Diameter Ø	Stärke Strength	Anwendungsbereich Field of application
M 943F /210 HP (fein, fine)	Ø 21 mm	0,16 mm	Keramik; ceramic
M 943F /250 HP (fein, fine)	Ø 25 mm	0,16 mm	Keramik; ceramic



Sinterdiamant-Scheiben / Superflex
Sinterdiamond-Discs / Superflex

Beschreibung: Sinterdiamant-Scheiben / Superflex

Description: Sinterdiamond discs / Superflex



Bestellnummer Order Number	Durchmesser Ø Diameter Ø	Stärke Strength	Anwendungsbereich Field of application
M 944F /080 HP (fein, fine)	Ø 0,8 mm	0,15 mm	Keramik; ceramic
M 944F /120 HP (fein, fine)	Ø 12 mm	0,15 mm	Keramik; ceramic
M 944F /160 HP (fein, fine)	Ø 16 mm	0,15 mm	Keramik; ceramic
M 944F /180 HP (fein, fine)	Ø 18 mm	0,15 mm	Keramik; ceramic
M 944F /210 HP (fein, fine)	Ø 21 mm	0,15 mm	Keramik; ceramic
M 944a /210 HP (mittel, medium)	Ø 21 mm	0,17 mm	Keramik; ceramic
M 944b /210 HP (mittel, medium)	Ø 21 mm	0,23 mm	Keramik; ceramic

Notiz / Notice: