



Trepane / Trephines

rostfrei / stainless



Aussen-Ø
External Ø
1,8 2,0 mm



223 RS

HP	L: 44,5 mm	008 010
----	------------	---------

rostfrei / stainless



Aussen-Ø
External Ø
2,8 3,3 3,7 mm



224 RS

HP	L: 44,5 mm	018 023 027
RA L	L: 26,0 mm	018 023 027

rostfrei / stainless



Aussen-Ø / External Ø
2,4 2,8 3,1 4,0 5,0 6,0 7,0 8,0 mm



225 RS

HP	L: 44,5 mm	014 018 021
RA	L: 22,0 mm	014 018 021 030 040 050 060 070

rostfrei / stainless



Aussen-Ø
External Ø
6,0 mm



227 RS

HP	L: 44,5 mm	050
----	------------	-----



Trepane / Trephines

rostfrei / stainless

229 RS

Aussen-Ø / Outside Ø mm

	3,00	3,50	4,00	4,50	5,00	5,50	6,00	
RA L	L: 26,0 mm	020	025	030	035	040	045	050

Länge Arbeitsteil / Length Working Part: 10,0mm

rostfrei / stainless

229 RS

Aussen-Ø / External Ø

	6,50	7,00	7,50	8,00	8,50	9,00	9,50 mm	
RA L	L: 26,0 mm	055	060	065	070	075	080	085

Länge Arbeitsteil / Length Working Part: 10,0mm

rostfrei / stainless

229 RS

Aussen-Ø / External Ø

	10,00	10,50	11,00	11,50	12,00	12,50	13,00 mm	
RA L	L: 26,0 mm	090	095	100	105	110	115	120

Länge Arbeitsteil / Length Working Part: 10,0mm

rostfrei / stainless

229 L RS

Aussen-Ø / External Ø

	4,00	4,50	5,00	5,50	6,00	6,50	7,00	7,50	8,00	8,50	9,00 mm	
RA L	L: 26,0 mm	030	035	040	045	050	055	060	065	070	075	080

Länge Arbeitsteil / Length Working Part: 14,0mm

rostfrei / stainless

229 XL RS

Aussen-Ø / External Ø

	4,00	4,50	5,00	5,50	6,00	6,50	7,00	7,50	8,00	8,50	9,00 mm	
RA L	L: 26,0 mm	030	035	040	045	050	055	060	065	070	075	080

Länge Arbeitsteil / Length Working Part: 18,0mm